




Clarity Trial Mode

This software is delivered in a Trial mode (without a user code). Please leave the user code field empty during the Clarity installation. The hardware key must be inserted in the computer.

The **TRIAL** mode can be recognized by the **TRIAL** inscription in the header of the main **Clarity** window.

TRIAL mode enables a user to work in the **Clarity** station with full functionality for the duration of 30 days or 100 runs of the **Clarity**, whatever occurs first.

After the trial mode period has expired, the station will switch to the **DEMO** mode, where data acquisition and some other functions are disabled. The station in the **DEMO** mode is recognized by the **DEMO – Trial Expired** inscription in the red label in the main **Clarity** window.




Hardware key replacement

In case you have received new hardware key and you already entered the User Code in the **Clarity (Lite)**, please delete the *Clarity.sno* file from the **Clarity (Lite)** installation directory and restart the application.

In order to get fully functional software (Full mode) you need to insert a user code. Please follow instructions below:

Switching Clarity from Trial to Full mode

1. Open the **Clarity** station.
2. Use the **Help – User Code...** menu command (as shown in the picture) to open the **User Code** dialog.
3. Fill in the correct User Code sent to you by the **DataApex** company or your local distributor.
4. Use the **OK** button to accept the user code.
5. The **User Code** dialog will close and in just only a few seconds, the station will switch to the full mode.



Switching Clarity LITE from Trial to Full mode

The procedure is the same as when upgrading the **Clarity** station, but as the **Clarity LITE** version does not contain the main **Clarity** window, the **Help – User Code...** menu command must be invoked from the **Clarity LITE** window. The rest of the steps (from step 2) is the same as in the previous procedure.

